

## Preliminary experiments in FYS-KJM4740

**Objective:** Elucidation of some basic NMR concepts by practical NMR

The present NMR tutorial (room V170 in the Department of Chemistry) is an offer to students who want to learn some of the NMR concepts by a practical approach. The experimental results will be available on the home page of the course during January. All students must write up a short report based on the experimental findings.

### Exercise 1

Experiments on water confined between glass beads:

Practical: Insert the sample in correct position

**Acquisition mode** → Sequence → Load → FID

Parameters: VT 30C, DW 0.5, RD 3s, NS 4, RG 10, SI 1K, P90 1.95. The  $t_d$  (rf pulse duration) is varied.

**Start experiment: GO**

- a) Set  $t_p$  (0.5, 1.0, 1.5, 2.00, 2.50, 3.0, 3.50, 4.00, 4.50 and 5.00  $\mu$ s) and determine the initial signal intensity  $I$  of the FID vs  $t_p$ . What is the strength  $B_1$  (gauss) of the rf-field? Discuss your results.
- b) Acquire an FID (real/imaginary) for different rf-frequency offsets ( $O1$ ). Plot the FID for different  $O1$  and discuss your results with respect to the rotating frame of reference.  
Set: RD 3s, VT 30C, DW 0.5, SI 1K, P90 1.95, NS 4, RG 10).  $O1$  is varied.
- c) Fit a single exponential function and a Gaussian function to the on-resonance FID with  $t_p = 1.95\mu$ s. Discuss your results.
- d) Acquire an FID for different repetition delays  $t_{RD}$  (= 3s, 2s, 1s and 0.5s). Plot the initial signal intensity  $I(0)$  of the FID as a function of  $t_{RD}$  and discuss your results.  
Set: DW 2, SI 1K, VT 30C, P90 1.95, NS 4, RG 10, on resonance, vary  $t_{RD}$ .

## Exercise 2

Inversion Recovery experiments on a series of CuSO<sub>4</sub>-solutions (Concentrations: 20mM, 10mM, 5mM, 2.5mM, 1.25mM and 0mM).

Practical: Insert the sample in correct position

**Acquisition mode** → Sequence → Load → INVREC

Set: **RD 1s, RG 1, VT 30, NS 4, DW 0.1s, SI 4K**

**Write: .T1** → **T1-4740** → open → OK → write a file name.

A T1-plot will appear on the screen when the experiment is finished.

Plot T1 vs CuSO<sub>4</sub>-concentration and discuss your results.

## Exercise 3

CPMG experiments on a series of CuSO<sub>4</sub>-solutions (Concentrations: 20mM, 10mM, 5mM, 2.5mM, 1.25mM and 0mM).

Practical: **Acquisition mode** → Sequence → Load → CPMGF

**Set : NS 4, RD 1s, RG 1.5, VT 30, P90 = 2.0 (μs), TAU 100 (μs), NECH 4K.**

**Start: GO.**

**When finished;** Write T2.

A plot will appear on the screen. Make a “screen dump” (Press “Print Scrn” and paste the image into “Paint”). You will use the value of T<sub>2</sub> in the report.

Plot T2 vs CuSO<sub>4</sub>-concentration and discuss your results.

**Write a report on your experimental work (include relevant Figures)!**

## Practical NMR

### A Simplified User Manual for the Maran Ultra NMR Instrument (Non-expert User)

#### Initiate experiment

1. Insert a sample in correct position
2. **Acquisition mode:** sequence → load → FID.EXE
3. Set RD to  $\approx 5T_1$ , i.e., **RD 1.0s, P90 = 2.15** ( $\mu\text{s}$ ), **VT 30** ( $^{\circ}\text{C}$ ) and **RG 1, DW 0.1, SI 4K**.

#### Optimizing relevant parameters

4. **Command** → *Auto OI* (wait)
5. **Command** → *Auto RG* (wait). Use this value of RG throughout if not otherwise stated in the text.
6. **Command** → *Auto P90* (wait)

#### T<sub>1</sub> measurement

7. Insert the sample in correct position
8. **Acquisition mode** → Sequence → Load → INVREC
9. Set: **RD 1s, RG 1, VT 30** and **NS 4, DW 0.1s, SI 4K**
10. **Write: .T1** → *T1-4740* → open → OK → write a file name.
11. A T<sub>1</sub>-plot will appear on the screen when the experiment is finished. Make a “screen dump” (Press “Print Scrn” and paste the image into “Paint”). You will use the value of T<sub>1</sub> in the report.

#### T<sub>2</sub>-measurement

12. Insert the sample in correct position
13. **Acquisition mode** → Sequence → Load → CPMGF
14. Set: **RD 1s, RG 1, VT 30** and **P90 = 2.15** ( $\mu\text{s}$ ), **tau 100** ( $\mu\text{s}$ ), **NECH 4K**
15. **Write: GO**
16. **Store File:**
17. Write T<sub>2</sub>. A plot will appear on the screen. Make a “screen dump” (Press “Print Scrn” and paste the image into “Paint”). You will use the value of T<sub>2</sub> in the report.

#### (Miscellaneous)

##### Calculations using the NMR program WinFit

18. Open WinFit: C:\ProgramFiles\Resonance\WinFit\nmrfit (double click)
19. Retrieve a file: File → Open → Select file → data → (Select Fit Options) → number of exponents/DC-offset (if possible) → Auto Initialize → Fit → Data (Read out values) or print, i.e.; File → Print → ....