

Laboratorieøvelse i FYS3180/4180

Eksperimentelle metoder i fysikk

Oslo Cyclotron Laboratory

<http://ocl.uio.no/>

Formål:

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
Fig. 3: Nal responsfunksjonen.

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AR		ARithmetic's (+-/*)
CA		Get info and give a0, a1 and a2 and dimension
CE		Change between Channel and Energy display
CL		Calibration Lines
CO		COmpress spectrum
CR		Call Cursor and get info
DS		DiSplay spectrum
DX	M1 M2	Low/high markers for Display of X axis
DY	M1 M2	Low/high markers for Display of Y axis
DZ	M1 M2	Low/high markers for Display of Z axis
EX		EXpand display using cursor (mouse)
FD		Fit Data points with polynomial, exponential, etc.
FO		FOLD spectrum with the response matrix
FT	N	N=1-15: set-up for N peaks and do fit
GR		Get Response matrix into working matrix 1 or 2
HE		List commands
IC		Select spectrum
LS		LiSt files in current directory (ls from UNIX)
DM		Display Matrix as spectra
OS		OVERlay of spectra
PO		POLynom fit to set of data
PF		Peak Find on spectrum or matrix
PM		Project Matrix down on X or Y axis
RE		REAd spectrum
RM		Make Response Matrix
SC	N	SCale for Y axis: N=1/2/3=lin/quadratic/log10
SD		Sirius Display, display acquisition spectra
ST		STop and exit Mama
SU	M1 M2	SUM counts between channel M1 and M2 using cursor
UN		UNfold spectrum with response matrix
WR		WRite spectrum with calibration and dimensions