Synthesis of BaTiO₃ by hydrothermal methods and via a heterometallic alkoxide.

Introduction:

During the last decades the use of functional oxides for new applications has led to increased demand regarding purity and homogeneity. Traditional methods used for synthesis of heterometallic oxides usually employ high temperature and will often result in inhomogeneity at a microscopic level. Other sol-gel methods, like the citrate route, do circumvent these problems, but may lead to contamination with carbon, often in the form of carbonates. This is especially important for elements, which easily form carbonate.

When utilizing the materials in microelectronic devices it is often necessary to synthesize the materials as thin films (e.g. 100-500 nm thickness for ferroelectrical capacitors). Methods used for producing films in this thickness range are e.g. sputtering, PLD (pulsed laser deposition) and CSD (chemical solution deposition). In the CSD method a solution (or sol) of metalorganic compounds with the same stoichiometry as the desired product is used. CSD is much used as the method does not require much equipment and there is no need for time consuming optimization of parameters in order to obtain the correct stoichiometry. For practical application, the crystallization temperature is important in order to avoid damage to other components during the heat treatment. The crystallization temperature is reduced when reducing the diffusion distance needed for the atoms to organize into a crystalline arrangement. Therefore, homogeneity on a nanometer scale is an advantage. This is not accomplished by traditional sol-gel methods where several precursors are used.

Multimetallc alkoxides are molecular species which contain two or more metal atoms. These molecules then already have a homogeneous distribution of the metal atoms on a nanometer scale, which is maintained during a decomposition process. Crystallization often takes place at temperatures from 450°C to 550°C. Thus, metal alkoxides are attractive precursors for metal oxide synthesis.

Ferroelectricity was observed for the first time in BaTiO₃ ($T_C=130$ °C). This material takes the perovskite structure and below T_C it is tetragonal with a dipolar moment along the c -axis. Above the Curie temperature the material is cubic with zero net dipolar moment. The material transforms to structures with lower symmetry below room temperature. In tetragonal BaTiO₃ the potential of Ti⁴⁺ along the c -axis form a double-well. An electrical field opposite the polarization direction may pull Ti⁴⁺ over the potential barrier, thereby reversing the polarization. BaTiO₃ is a technologically

important dielectric material and is used e.g. in capacitors ($\epsilon_r = 2000-10000$) and PTC-resistors (positive temperature coefficient). All ferroelectrics are also piezoelectric, but other piezoelectric materials are more important than BaTiO_3 .

Particles formed by crystallization of alkoxide-gels are very small. This makes the phase transition at 130°C by cooling energetically unfavourable. The powder resulting from the synthesis is therefore in a metastable cubic modification and must be heated to $800-900^\circ\text{C}$ in order to obtain traces of the tetragonal phase forming by cooling. In order to achieve 100% conversion, heating to above 1000°C is necessary.

As it is the tetragonal phase which is ferroelectric and therefore the technologically important material, the possibility of synthesizing tetragonal BaTiO_3 directly at low temperature has been pursued. The challenge is to obtain sufficient crystal growth at low temperature in order to avoid the tetragonal-cubic phase at 130°C . This may be possible by using hydrothermal synthesis, where it is possible to crystallize metal oxides from 150°C .


In this exercise we will use X-ray powder diffraction to investigate structure (symmetry) and phase purity.

Synthesis:

Alkoxide:

We want the organic content to be as low as possible in order to avoid carbon contamination. Small alkoxide groups like methoxy and ethoxy result in either unstable alkoxides or alkoxides with low solubility. Iso-propoxy is a much used alkoxy group, and we are going to use this for our reactions. A heterometallic alkoxide may be prepared by reacting the homometallic alkoxides under reflux. Titanium isopropoxide is readily available and therefore we will synthesize only barium isopropoxide. This is an oligomeric compound, $1/m[\text{Ba}(\text{OPr}^i)_2]_m$, which is dissolved during the reaction with Ti-isopropoxide. This is also a good indication that a larger heterometallic species is formed, as these generally have higher solubility.


Alkoxides of electropositive metals may be synthesized by direct reaction between metal and the alcohol (Lewis base), e.g.:


The reaction rate increases with the acidity of the alcohol, decreasing steric hindrance, basicity of the metal and the solubility of the product. All alkali metals react very fast with light alcohols to give pure well defined products. The alkaline earth metals give

more complex products and the least electropositive metals need a catalyst in order to speed up the reaction rate (e.g. Mg). Barium is, however, very electropositive and needs no catalyst. When a very soluble product is formed cooling is necessary. Barium isopropoxide is, however, insoluble and forms a barrier between metal and alcohol which slows down the reaction. A suitable reaction rate between barium and isopropanol at room temperature is obtained by using fairly small particles, 0.2 – 0.5mm. Note that barium is generally very reactive toward hydroxyl containing compounds, such as water. **All handling of barium must therefore be done in inert atmosphere.**


Synthesis of alkoxides of metals which do not react directly with alcohols may be performed by alcolysis of metal amides. The starting compound is here a metal halide:

LiCl is removed by filtering and the alcolysis is readily done by adding alcohol to the amide:

Another method is to react metal halides with alkoxides of alkali metals, e.g.:

Reactions often give less than 100% yield and purification e.g. by sublimation or recrystallization must be used to enrich the alkoxide.

Things to consider during the exercise:

You must now synthesize and handle compounds which are very air sensitive. Directions of how to use the vacuum line will be given. However, not all steps may be covered; it may not be specifically noted that a flask needs to be evacuated on the vacuum line before connecting it to the rest of the system. **It is therefore important to constantly consider whether air may enter the system in any way and get into contact with other compounds on the vacuum line.** You must therefore develop routines and sequences for opening and closing of valves for various operations e.g. putting an empty flask on the vacuum line. Remember that air diffuses from one place to another. **It is therefore important to close the access to other flasks containing air sensitive compounds before evacuating a flask containing air. It is also important to avoid access to flasks containing solvents when you evacuate the**

line. Otherwise the content may boil and enter the line. All ground glass joints (NO: slip) must be secure by clamps or elastic bands.

Readying the vacuum line:

In order to protect the vacuum pump there must always be liquid nitrogen in the cold trap before the pump. Start the pump with valves 1 and 2 closed (refer to drawing). If there is a gas ballast valve this should be left open until the pump is warm. This increases the life time of the pump. Open valves 2, 4, 8 and 10 in order to evacuate the system. After a couple of minutes valve 2 may be closed, then nitrogen is released into the system from the cylinder (until bubbles are observed in the mercury container)???. Then evacuation is repeated.


Figure 1: Illustration of the vacuum line for synthesis under inert conditions.

Solvent drying:

When synthesizing compounds which are sensitive to Lewis bases it is important to dry all solvents used. Liquids as hexane, benzene, THF, di-ethyl ether etc. do not react with sodium metal. These may therefore be pre-dried by refluxing with 5x5mm pieces of sodium for 1 hour. In order to remove the last traces of water and dissolved gasses, 10g benzophenone per liter is added and the mixture is refluxed for until a blue/violet

solution is obtained. The colour indicates the presence of a radical and the solution is therefore free of reactive species. The solvent may then be distilled.

Alcohols, on the other hand, react with sodium and must be dried by refluxing with e.g. CaH_2 or a corresponding metal alkoxide. After 12 hours the solvent may be distilled.

Solvents which do not react with sodium metal are stored over 5x5mm sodium bits or wire. All solvents are stored under nitrogen in a schlenk-flask. Nitrogen flow is used when taking solvent by a pipette in order to avoid uptake of moisture from the air.

For this exercise you will be receive dry solvents.

Synthesis of $1/m[\text{Ba}(\text{O}^i\text{Pr})_2]_m$

Before starting the laboratory responsible will show you how to operate the glove box.

Please remember always to close the air lock from the inside when you work. This is easy to forget if e.g. more equipment has to be brought into the glove box.

Take a clean, dry round-bottomed flask, a funnel with ground glass joint and a glass stopper with some vacuum grease and a spatula into the glove box. Ba(s) is found inside the glove box. Weigh 0.5g Ba, put on the stopper, secure with a string and take the equipment out of the glove box. Evacuate a T-piece and place the flask in this.

An important thing to remember when working with liquid/solid reactions which are fast and/or releases heat is to always dissolve or disperse the solid in an inert solvent. If you add a reactive liquid directly to a solid it will produce a large amount of heat just where the liquid is in contact with the solid without any possibility of transporting the heat away. Decomposition of the product may occur. Therefore, 20 ml dry hexane and a bar magnet is added to the flask containing the barium metal. Calculate the amount of isopropanol needed and add a surplus to the flask under stirring. After a short while, evolution of hydrogen will start and then $1/m[\text{Ba}(\text{OPr}^i)_2]_m$ will precipitate. Continue stirring until all barium has reacted.

Synthesis of $[\text{BaTi}(\text{OH})(\text{O}^i\text{Pr})_5(\text{Pr}^i\text{OH})_3]_2$:

Synthesis of heterometallic alkoxides are often performed by simply reacting the monometallic alkoxides¹. As this is a Lewis acid-base reaction, reactivity decreases with decreasing difference in Lewis acid/base strength of the reactants.

Task:

Write the reaction equation for formation of $[\text{BaTi}(\text{OH})(\text{O}^i\text{Pr})_5(\text{Pr}^i\text{OH})_3]_2$. What must be present in addition to $1/m[\text{Ba}(\text{OPr}^i)_2]_m$ og $\text{Ti}(\text{O}^i\text{Pr})_4$?

We chose to weigh the amount of $\text{Ti}(\text{O}^i\text{Pr})_4$ instead of using a pipette as this will be more precise. This must be done in the glove box. Calculate how much $\text{Ti}(\text{O}^i\text{Pr})_4$ is needed in order to obtain a ratio between Ba and Ti of **1:0.98**. This is done in order to ensure a perfect stoichiometry in the final product. *Can you think of why we may obtain stoichiometry by using non-stoichiometry here?*

In the glove box you will need:

- A couple of Pasteur pipettes with rubber ball.
- A flask containing **dry** 2-propanol, secured by strings.
- A weighing container
- The flask containing $1/m[\text{Ba}(\text{OPr}^i)_2]_m$ secured with strings.

Weigh the calculated amount of $\text{Ti}(\text{O}^i\text{Pr})_4$ into the weighing container and add this to the flask. Flush with 2-propanol. Be careful not to spill anything and rinse well. Replace the stopper, take the equipment out of the glove box and mount the flask on the cooler, which previously has been mounted on the vacuum line via a T-piece, and evacuated. Reflux for three hours. Note the dissolution of barium isopropoxide after a short while.


Figure 2: $[\text{BaTi}(\text{OH})(\text{O}^i\text{Pr})_5(^i\text{PrOH})_3]_2$

As we have used a small surplus of $\text{Ti}(\text{O}^i\text{Pr})_4$ you will notice that the solution does not become entirely clear. As the molecule produced does have an exact 1:1 stoichiometry between Ba and Ti we may now ensure perfect stoichiometry by filtering off the surplus of $1/m[\text{Ba}(\text{OPr}^i)_2]_m$.


When doing this you must take the flask off the line and into the glove box. Remember the strings! For filtering a 0.2µm syringe filter is used. You may leave the flask for precipitation for a few minutes in order to avoid clogging the filter. Place the piece of PVC hose on the syringe and suck up the content from the top, avoiding too much precipitate. Then switch to using the filter, and press the liquid through. This may be a bit hard. **Be careful that the filter is not pressed off or destroyed.** If the filter clogs up, it may help to pull the syringe back a bit.

When you have done this you have obtained a solution of $[\text{BaTi}(\text{OH})(\text{O}^i\text{Pr})_5(^i\text{PrOH})_3]_2$ -molecules in Pr^iOH and hexane.

Hydrolysis:

In order to transform the solution to a gel-like state the alkoxide must be hydrolysed. During hydrolysis a three dimensional framework of oxygen bridges is formed.

Alkoxides hydrolysis in contact with water:

At the same time a condensation process occur:

The degree of hydrolysis is defined as:

$$h = [\text{H}_2\text{O}]/[\text{M}(\text{OR})_x] \quad \text{Eq. 7}$$

and depends on the purpose of the produced material. The transition from a sol to a gel is a very complex process and may involve many more steps than described here. By controlling

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