Germaniumdetektorer og spektroskopianalyse

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English translation

Learning goal: You will learn about semiconductor detectors, in particular Ge-detectors. Such detectors are used for precise determination of the intensity and energy of γ -quanta, commonly referred to as γ spectroscopy.

In the exercise we will not perform the actual analysis of γ spectra, but we will look into how the results can be used to identify an unknown sample, and how we can determine its absolute amount.

The theory behind this exercise you can read about here:

- ! Web-page about γ radiation, written by Liv Stavsetra (in Norwegian): http://www.uio.no/studier/emner/matnat/kjemi/KJM5900/h04/Ekstramateriale/Gammastraling_v02_Sep2004.html
- ! Chapter 7 in the course book (page 95-125) gives a general description of how to measure radioactivity. Information particularly about γ measurements you will find in the sub chapters 7.4 Semiconductor Detectors (page 106), 7.5 Choice of Detectors (page 110) and 7.6 Spectrometry (page 113).
- ! Chapter 6.4 " γ Radiation" (page 85) in the course book describes how γ radiation interacts with matter.

Useful tables and references:

- ! Table with standard γ sources is on page 15 in the booklet to the Karlsruher Nuklidkarte.
- ! "Gammas" database, which is installed on the PCs in the course laboratory.
- ! Table of Isotopes (LBNL/Lund): <http://ie.lbl.gov/education/isotopes.htm>.

Exercise purpose

In this exercise you are to do an efficiency calibration of a Ge-detector. The efficiency to any detector is dependent on quantities such as the γ -rays energy, the geometrical shape and size of the source, and the distance between the source and detector.

In order to make an efficiency calibration we use sources with known amounts of activities. The sources are measured and then we can calculate the counting efficiency: If the source activity is A , the number of counts obtained with a given detector is R , then the counting efficiency is given by:

$$\varepsilon = \frac{R}{A}$$

To calibrate a γ detector we must also take into consideration that a γ ray will not necessarily be emitted every time the source disintegrate. The frequency γ rays are emitted can be found in tables and is given as the γ intensity, usually symbolised with I_γ . The formula we use to calculate the efficiency then looks like this:

$$\epsilon_E = \frac{R_E}{A \cdot I_\gamma}$$

Notice that we have indexed the ϵ and R with an E to indicate that these are energy dependant. Since all the absorption processes for γ radiation is strongly dependent on the γ -ray energy, the detector efficiency will also be dependent on energy. Thus, to calibrate a γ detector we must use a range of sources with energies that cover the range of energy we are planing to measure in. The efficiency calibration will then look like something in Figure 1. The most common range to measure γ radiation in is from about 50 keV and up to about 1-2 MeV. From Figure 1 we can se that the range between about 200 keV and up to 2 MeV is nearly linear in double-logarithmic scale. Below 200-250 keV this is not true.

In this exercise we have measured a series of standard γ -ray sources with a HPGe-detector (the HP in front of the element symbol Ge is short for *high purity*). The number of counts measured you shall use to plot an efficiency calibration curve for the detector. After that you shall use the efficiency curve to determine the amount of an unknown metal in an sample which has been irradiation with neutrons.

Exercise 3A: The Ge-detector

Read in the course book about Ge-detectors. If you want to, you can also try to find information on the net or in other course books (some you'll find in the course laboratory). Write a short summary, about half a page, about what a Ge-detector is, it's properties and how it works.

Exercise 3B: Efficiency Calibration Curve

Make an efficiency calibration curve, where the counting efficiency is given as a function of γ energy. In table 1 the measured nuclei is listed with γ energy, number of counts, uncertainty in the number of counts and the counting time. The measurements was performed 27th March 2003.

A table of the γ intensities, I_γ , you can find at page 15 in the Karlsruhe Nuklidkarte booklet.

To calculate the efficiency you need to know the amount of activity in the sources. Table 2 lists the amount of activity at the time of purchase. You must correct these number based upon the age and half life of the sources.

Exercise 3C: Unknown sample

In table 3 the analysis results of the unknown, n-irradiated sample is shown. Use the results to :

- ! Determine which metal the sample consist of. Use the identified γ -ray energies to do this. Databases is available in the counting room (contact your supervisor if you want a copy to install on your home computer).
- ! Determine the absolute amount of the metal. Use the counts measured in table 3 and the efficiency curve you made in exercise 3B. Use at least two different γ lines to do this, we recommend $E_\gamma = 222.08$ keV and $E_\gamma > 1$ MeV. The neytron flux was 10^5 n/s cm² during the irradiation, and the irradiation time was 18 days. The sample was measured directly after the irradiation.

Your report

In this exercise, the report is not scrutinised as carefully as in exercise 2. But make sure that your report is neat and tidy. Include the following in your report:

- ! Describe what a Ge-detector is and how it works.
- ! Tabulate values and intermediate results you used to plot the efficiency curve. Remember to calculate the uncertainty in the individual points!
- ! You are free to choose if you want to plot your data by use a computer program (e.g. Origin or MS Excel) or draw them manually.
- ! The determination of the unknown, n-activated metall and its amount.