LAB FRM2041

29. MAI 2006

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STREPTOCOCCAL GROUPING LATEX KIT (PRO-LAB)

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**PROLEX™ STREPTOCOCCAL
GROUPING LATEX KIT**
(for *in vitro* diagnostic use)

INTENDED USE

Prolex™ streptococcal grouping latex kit provides a rapid method for the serological identification of groups A, B, C, D, F and G of the Lancefield groups of streptococci growing on agar plates.

SUMMARY AND EXPLANATION

Clinical, epidemiological and microbiological studies have conclusively shown that the diagnosis of streptococcal infections based on clinical symptoms always requires microbiological confirmation (4). Beta-haemolytic streptococci are the most frequently isolated human pathogens among the representatives of the genus *Streptococcus*. Nearly all the beta-haemolytic streptococci possess specific capsular carbohydrate antigens (streptococcal group antigens). Lancefield showed that these antigens can be extracted in soluble form and identified by precipitation reactions with homologous antisera. Different procedures for extraction of streptococcal antigens are currently in use (1,2,6,7,10,11). The Prolex™ streptococcal grouping kit is based on the liberation of specific antigen from bacteria cell wall by modified nitrous acid extraction. The extracted antigen in conjunction with latex agglutination offers a rapid, sensitive and specific method for identification of streptococcal groups A, B, C, D, F and G from primary culture plates.

PRINCIPLE OF THE TEST

The Prolex™ streptococcal grouping method involves chemical extraction of group specific carbohydrate antigens using specially developed nitrous acid extraction reagents. The extraction reagents 1 and 2 provided in the kit contain a chemical substance able to extract the streptococcal group specific antigens at room temperature. Extraction 3 contains a neutralizing solution. The neutralized extracts can be easily identified using latex particles sensitized with purified group specific rabbit immunoglobulins. These latex particles agglutinate strongly in the presence of homologous antigen and will not agglutinate when homologous antigen is absent.

REAGENTS

Each kit is sufficient for 60 streptococcal grouping tests. Materials are supplied ready for use.

Latex Suspensions: Six vials each containing 3.0 ml of latex particles coated with purified rabbit antibodies to Group A, Group B, Group C, Group D, Group F or Group G streptococci. The latex particles are suspended in phosphate buffer pH 7.4 containing 0.1% sodium azide as preservative.

Polyvalent Positive Control: One vial containing 2 ml of ready to use polyvalent antigens extracted from inactivated streptococci of Lancefield groups A, B, C, D, F and G. The strains for antigen preparations are ATCC strains listed in the section "MATERIALS REQUIRED BUT NOT PROVIDED".

Extraction Reagent 1: One dropper bottle containing 3.2 ml of extraction reagent 1 with 0.1% sodium azide as preservative.

Extraction Reagent 2: One dropper bottle containing 3.2 ml of extraction reagent 2.

Extraction Reagent 3: Two dropper bottles each containing 8 ml of extraction reagent 3 with 0.1% sodium azide as preservative.

Canada • 905-731-0300 Fax 905-731-0206
20 Mural Street, Unit #4, Richmond Hill, Ont. L4B 1K3

All components of the kit are available separately:

- Latex Suspension Group A PL.051
- Latex Suspension Group B PL.052
- Latex Suspension Group C PL.053
- Latex Suspension Group D PL.054
- Latex Suspension Group F PL.055
- Latex Suspension Group G PL.056
- Extraction Reagent 1 PL.057
- Extraction Reagent 2 PL.058
- Extraction Reagent 3 PL.059
- Polyvalent Positive Control PL.060

PRECAUTIONS

1. Do not use reagents after expiry date shown on product label.
2. Some reagents contain sodium azide. Sodium azide can react explosively with copper or lead if allowed to accumulate. Although the amount of sodium azide in the reagents is minimal, large quantities of water should be used when flushing used reagents down the sink.
3. Extraction reagents contain caustic agent, wash with soap and copious amounts of water immediately in case of skin contact, and flush generously (at least 15 minutes) with water in case of eye contact.
4. Safety precautions should be taken in handling, processing and discarding all clinical specimens as a pathogenic organism may be present.
5. The kit is intended for *in vitro* diagnostic use only.
6. The procedures, storage conditions, precautions and limitations specified in these directions must be adhered to in order to obtain valid test results.

STABILITY AND STORAGE

All kit components should be stored at 2-8°C. Do not freeze. Reagents stored under these conditions will be stable until the expiry date shown on product label.

SPECIMEN COLLECTION AND PREPARATION OF CULTURES

For specific procedures regarding specimen collection and preparation of primary cultures refer to a standard microbiology textbook. In general, a fresh (18-24 hr.) Gram positive beta-haemolytic (5% sheep blood agar) isolate of streptococcal colonies is assumed. Four large colonies should be adequate for grouping; however, if the colonies are minute, an increased number of colonies (100pfu) should be used.

MATERIALS SUPPLIED

Latex suspensions for each of Group A, Group B, Group C, Group D, Group F and Group G streptococci.
Polyvalent positive control containing polyvalent extract representing antigens from streptococcal groups A, B, C, D, F and G.
Extraction Reagents 1, 2 and 3.
Disposable cards with eight test circles labelled 1 through 8.
Disposable mixing sticks.

PRODUCT CODE PL.050 (60 tests)

MATERIALS REQUIRED BUT NOT PROVIDED
Inoculating loops, Pasteur pipettes, Borosilicate glass test tubes 12mm x 75mm, Timer, Vortex mixer (optional).

The following ATCC strains are recommended for use in quality control.
Streptococcus pyogenes group A (ATCC# 19615),
Streptococcus sp. group B (ATCC# 12386),
Streptococcus sp. group C (ATCC# 12388),
Enterococcus faecalis group D (ATCC# 19433),
Streptococcus sp. Type 2, group F (ATCC# 12392),
Streptococcus sp. group G (ATCC# 12394)

TEST PROTOCOL

- All components should be at room temperature (22-28 °C) prior to use.
1. Label one test tube for each specimen.
 2. Add 1 drop of Extraction Reagent 1 to each tube.
 3. Select 4 large beta-haemolytic colonies with disposable loop and suspend them in the Extraction Reagent 1. If colonies are minute, pick several well isolated colonies to be tested such that Reagent 1 solution becomes turbid. In all cases the streptococcal colonies should be picked from an area which contains the least amount of contamination.
 4. Add 1 drop of Extraction Reagent 2 to each tube.
 5. Mix the reaction by vortexing the tube or by tapping the tube with a finger for 10-15 seconds. Incubate at room temperature for two minutes.
 6. Add 5 drops of Extraction Reagent 3 to each tube. Mix the reaction as in step 5.
 7. Dispense one drop of each latex suspension onto separate circles on the test card.
 8. Using a pasteur pipette, place one drop of extract beside each drop of latex suspension.
 9. Mix the latex and the extract with sticks provided using the complete area of the circle. A new stick should be used for each reagent.
 10. Gently rock the card allowing the mixture to flow slowly over the entire test ring area.
 11. At two minutes, under normal lighting conditions, observe for agglutination.

QUALITY CONTROL PROCEDURES

Routine quality control procedures for each Prolex™ lot involve testing of the kit components (latex suspensions, polyvalent positive control and extraction reagents) with extracts of each streptococci groups A, B, C, D, F and G using ATCC strains listed in the section "MATERIALS REQUIRED BUT NOT PROVIDED". In addition, each latex suspension is tested for absence of cross-reactions against extracts of the following organisms: *Escherichia coli* (ATCC #25922), *Klebsiella pneumoniae* (ATCC #13883), *Staphylococcus aureus* (ATCC #25923), *Hemophilus influenzae* type b (ATCC #10211) and *Streptococcus pneumoniae* (ATCC #49619). However, the following procedures are recommended to check the performance of the reagents:

U.S.A. # 800 522 7740 Fax 800 332 0450
Austin, Texas, 78758

U.K. # 0151 353 1613 Fax 0151 353 1614
Neston, Wirral, L64 3UJ

1. The positive control is used to check the performance of the individual latex reagents. The latex reagent should show obvious agglutination with the positive control. The positive control is not used to demonstrate the specificity of the latex reagents nor to ensure that the extraction step was performed correctly and is functioning.
2. The extract from known strain should agglutinate with homogenous latex reagents. Refer to the list of recommended ATCC reference strains to be used in the "MATERIALS REQUIRED BUT NOT PROVIDED" section.
3. As test of absence of autoagglutination the latex reagents should not show agglutination with normal saline solution.

INTERPRETATION OF RESULTS

Usable results: A significantly rapid strong clumping of the latex particles to form an agglutination pattern in only one of the latex reagents indicates specific identification of the streptococcal isolate. Weak or/and granular reaction with more than one latex reagent should be ignored. A weak reaction with single latex reagent should be repeated using a heavier inoculum. The repeated test is considered positive if a visible agglutination occurs with only one of the latex reagents. Figure 1 illustrates a suggested scheme for grouping streptococci.

Nonusable results: A milky appearance without any visible agglutination of the latex particles.

EXAMINATION OF THE PROCEDURE

1. False negative or false positive results can occur if inadequate amounts of culture or extraction reagents are used.
2. The kit is intended for use in identification of beta-haemolytic streptococci. If an alpha or non-haemolytic streptococci is identified, the identification should be confirmed by biochemical tests (5,9). (Refer to suggested scheme for grouping streptococci).
3. False positive reactions have been known to occur with organisms from unrelated genera, e.g. *Escherichia coli*, *Klebsiella* or *Penicillium* (3,8). These are likely to non-specifically agglutinate all latex reagents.
4. Some strains of Group D streptococci have been found to cross react with Group G antisera, this strain may be confirmed as Group D by the bile-esculin test.
5. Enterococci can be differentiated from Group D streptococci by biochemical tests.

EXPECTED VALUES:

Among the 578 isolates tested by the Prolexin at two clinical trial sites they were 422 Group A, 133 Group B, 43 Group C, 32 Group D, 19 Group F, 119 Group G and 90 non-groupable strains. These specimens were from heterogeneous hospital patients and outpatients and involved isolates from throat swabs, vaginal swabs, cervical swabs, urethral swabs, nose swabs, penile swabs, hip incision swabs, pus samples, respiratory samples, blood cultures, and urine. No attempt was made to relate the results to age, sex or any other epidemiological factors assuming that these isolates represent a sampling of any group of patients.

PERFORMANCE CHARACTERISTICS

A. Cross-reactivity studies:

Prolexin kit was tested for cross-reactivity using 33 ATCC reference strains. The kit successfully grouped all streptococci containing Lancefield groups A, B, C, D, F and G (N=16). No cross-reactivity was observed during the testing of other streptococcal strains (n=7); nor of other non-streptococcal organisms (n=10).

B. Clinical performance studies:

Prolexin kit performance was evaluated at a microbiological centre in Oxford, England (data on file at Pro-Lab Inc., Richmond Hill, Ontario, Canada). In this study 468 primary cultures were tested by the Prolexin kit and an alternative grouping kit. Overall with 457 of 468 isolates tested (97.6%). Anomalous results (n=16; Thirteen of the 16) anomalous results agreed after retest which included 1 group A, 2 group B, 1 group D, 1 group F, 5 group G and 1 non-groupable strain. Two of the 3 disagreed isolates were further identified as non-beta haemolytic strains. The third disagreed isolate was grouped as group D with the Prolexin kit and as group F by the alternative kit. This strain gave a positive streptococcal result with the alternative kit following alternative grouping kit after retest of anomalous results occurred with 463 of 468 isolates tested (99.1%). The 468 isolates used in this study included 127 group A, 52 group B, 30 group C, 28 group D, 9 group F, 107 group G and 75 non-groupable strains.


2. A second performance study was carried out at a Health Centre in Ontario, Canada. In this study, 111 primary cultures were included (110 tested, 1 untested). All the strains were originally grouped by Lancefield precipitation reactions. All group D were further biochemically confirmed using a BE (bile esculin) and PYR (pyrazinole) and/or discase assay protocol. The primary cultures were tested in parallel using the Prolexin Streptococcal grouping kit and an alternative grouping kit. In this study, the overall agreement between Prolexin and Lancefield results occurred with 109 of 110 isolates tested (99%), while overall agreement between the alternative kit and Lancefield results occurred with 106 of 110 isolates tested (96.3%). The 110 primary isolates used in this study included 15 group A, 40 group B, 13 group C, 4 group D, 11 group F, 12 group G and 15 non-groupable strains.

REFERENCES

1. Ederer, G.M., Hermann, M.M., Bruce, R., Maisten, J.M. and Chapman, S.S. (1972). Rapid Extraction Method with Tronase B for Grouping Beta-Haemolytic Streptococci. *Appl. Microbiol.*, 23, 285.
2. El-Kholi, A., Wannamaker, L.W. and Krause, R.M. (1974). Simplified Extraction Procedure for Serological Grouping of Beta-Haemolytic Streptococci. *Appl. Microbiol.*, 28, 836.
3. Elliot, S.D. and Tait, J.Y. (1978). The Type-Specific Polysaccharides of Streptococcus suis. *J. Exp. Med.*, 148, 1699.

4. Facklam, R.R. (1980). Streptococci and Aerococci, Ch. 8 in *Manual of Clinical Microbiology*, 3rd Ed., Edited by Tenover, E.H., Balows, A., Hansten, W.J., and Tarran, J.P. American Society for Microbiology, Washington, D.C. page 88-110.
5. Facklam R.R. (1977). *Physiological Differentiation of Vitamins Streptococci*. *J. Clin. Microbiol.*, 5, 184.
6. Fuller, A.T. (1938). The Fermentable Method for the Extraction of Polysaccharides from Haemolytic Streptococci. *Brit. J. Exp. Path.*, 19, 130.
7. Maxted, W.R. (1948). Preparation of Streptococcal Extracts for Lancefield Grouping. *Lancet*, ii, 255.
8. Nowlan, S.S. and Deibel, R.H. (1967). Group O Streptococci I. Ecology, Serology, Physiology and Relationships to Established Enterococci. *J. Bact.*, 94, 291.
9. Peltz, D.N. (1984). Early Detection of Streptococci in Swabs by Latex Agglutination Before Culture. *J. Clin. Microbiol.*, 19, 432.
10. Kautz, L.A. and Randall, E. (1955). Use of Autochained Extracts of Haemolytic Streptococci for Serological Grouping. *Stanford Med. Bull.*, 13, 290.
11. Watson, B.K., Moellering, R.C. and Kunz, L.J. (1975). Identification of Streptococci: Use of Lysostyrene and Streptomyces albus Filtrate in the Preparation of Extracts of Lancefield Grouping. *J. Clin. Microbiol.*, 1, 274.

Figure 1 SUGGESTED SCHEME FOR GROUPING STREPTOCOCCI

*Some strains of group D have been found to cross react with group G antisera Harvey, C.L. and McLinnery, M.R. (1984) *Can. J. Clin. Microbiol.* 10(41)